

Invacare® | **Birdie® Family**

UPDATED



Birdie and **Birdie Plus**
170 kg now upgraded to
180 kg.

Birdie Compact 150 kg



Birdie

Birdie Compact

Birdie Plus

Transportability and patient comfort go hand-in-hand

With maximum space for the user, both the *Invacare Birdie*, *Birdie Plus* and the *Birdie Compact* offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. The lifters are designed to ensure that folding and unfolding can be carried out easily and without the need for tools. In addition, the lifter can be dismantled into two parts if needed – also without tools.



Yes, you can.®



Total user comfort

Maximizes space for the user even when the hoist is in the highest position. This gives a very comfortable lift and transfer since the user can easily be rotated in the hoist. Consequently, it makes the placing of the user into the bed/wheelchair much easier for the carer. The optimization of space creates a high level of comfort since the user's knees will not be close to the actuator of the hoist. You have, then, a very low risk of injuring the knees of the user.

Ease of handling: simple to fold, unfold and dismantle without tools

Handling the lifter has been made as smooth and easy as possible as the hoists can easily be folded and unfolded without the need for tools. Hoists take up minimal space when stored and are easy to push and transport. If required, the mast and boom can be easily detached from the base (no tools required).

TÜV-certification for maximum security

For total security, the hoists come with a snap-lock system, thus minimizing the risk of accident. The lifters are CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure

the hoists comply with EN ISO 10535 (the European standard for hoists) the hoist has been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

Features and Options



Long boom – high user comfort

The user can easily be rotated and positioned. No risk of injuries to the knees.



Easy to fold and dismantle

Ease of transport and storage. No tools required when folding / unfolding or when dismantling the hoist.



Snap lock

Optimal security to avoid risk of accidents.



Emergency stop, emergency lowering and integrated charger

Manual or electric emergency lowering (compact version features only electric lowering). To charge the hoist, simply plug in.

Features and Options



Only available on Birdie Plus:

- Even Easier manoeuvring - 125mm rear wheels
- Electrical emergency lifting
- Showing battery status
- Soft start and stop
- Covered electrical legspread - noise reduction



Click on spreader bar
Easy to change the spreader bar.



Mounted spreader bar
Correctly and safe positioning of the spreader bar.



Manual leg spread
Simply step on the pedal to adjust the leg spread.



Handle for leg spread
Handle can easily be attached to operate the manual leg spread.



Electrical leg spread
Available for Invacare **Birdie** Standard version.



2 point quick release spreader bars
Easy to attach and detach. Available widths: 35 cm, 45 cm, 55 cm.



Push bar and hand control
Superior ergonomics for carers.



Castors
100 mm or 75 mm castors available. (braking castors are shown)













Scale
If needed the scale is available approved, according to the medical standard EN 45501, certified and sealed in order to be used to medicate.



4 point quick release spreader bars
Easy to attach and detach. Available widths: 45 cm and 55 cm.

Technical data

					
Birdie	1240 cm	64 cm	42 kg	44,5 - 171 cm	10 cm
Birdie Compact	109 cm	52 cm	36 kg	52,5 - 161,5 cm	10 cm
Birdie Plus	1240 cm	64 cm	47 kg	44,5 - 171 cm	10 cm

					
Birdie	101 cm	180 kg	40 cycles	2,9 Ah	140 cm
Birdie Compact	84,5 cm	150 kg	40 cycles	2,9 Ah	107 cm
Birdie Plus	91 cm	180 kg	40 cycles	2,9 Ah	140 cm

All measures are taken with 75 mm castors and 2 point spreaderbar.
For 100 mm castors add 1,5 cm for height measures and 2 cm for width.

Frame colour



Light grey /
Dark grey

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2014 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Birdie Family - EU - 10/2014

